


**ResMed**

# Hospital Trolley


ResMed's Hospital Trolley is well suited to a variety of hospital environments, providing optimised convenience and uncompromised mobility.



ResMed's Hospital Trolley has been intelligently designed to transport the Astral™ 100 and Astral 150 life support ventilators around any hospital or clinical setting.

It features the option of a medical grade power strip, located in the centre console, which lets you connect multiple electrical devices and power them with just one cord. This eliminates unnecessary wires, providing greater convenience and ease-of-use when moving ventilated patients.

The trolley also has the capacity to hold compatible accessories, including an oxygen cylinder, an external humidifier, and a circuit support arm—offering you an all-in-one, comprehensive solution.



## Key features

- Designed to hold the device, plus its full range of compatible accessories.
- Sleek, stable and highly mobile, the trolley can be transported with ease. The wheels have been designed to glide over uneven surfaces, and they are lockable, so you can safely secure the trolley.
- The humidifier bracket makes it easy to hook on an external humidifier in one simple step.
- You can store any necessary extras in the optional storage bin for even more efficiency.

## Technical specifications

|                        |  |
|------------------------|--|
| <b>Height</b>          | 1038 mm (1657 mm with IV Pole)<br>(40.9") (65.2" with IV Pole) |
| <b>Base diameter</b>   | 511 mm x 504 mm (20.12" x 19.8")                               |
| <b>Shipping weight</b> | 31.5 kg (69.5 lb)  |
| <b>Load limit</b>      | 25 kg (55.12 lb)   |

## Product codes

|  |              |
|--|--------------|
| Hospital Trolley   | <b>27914</b> |
| Hospital Trolley Base Plate (Astral)<br>(included with Hospital Trolley) | <b>27928</b> |
| Power Strip: US  | <b>27931</b> |
| Power Strip: EU  | <b>27932</b> |
| Power Strip: ANZ   | <b>27933</b> |
| Power Strip: UK  | <b>27934</b> |

Discover more about ResMed  
at [ResMed.com](http://ResMed.com)

|   |              |
|---|--------------|
| Hospital Trolley Cylinder Mount<br>(optional) | <b>27930</b> |
| Hospital Trolley Storage Bin (optional)       | <b>27929</b> |
| Circuit Support Arm                           | <b>24917</b> |